



# FILLING NOZZLE

## – GAS PURGING AROUND FILLING NOZZLE

- For gas purging around filling nozzle for liquid
- Protects liquid to be filled from contact with surrounding air
- Can be used together with filling nozzle 30-030-xxx
- Fits nozzle holders for FMB210, FPC50, FP50 and FF20



SPECIFICATIONS			
Materials	AISI 316L + Viton		
	C (ID)	D (OD)	Part number
Dimensions in mm	Ø4,5	Ø 5,3	30-031-050
	Ø8,5	Ø10,1	30-031-090
O-ring for gas filling nozzle			54-020-116
Nozzle holder incl. arm for filling stand			30-700-053
Declaration of compliance - Type 2.1 acc.: DS/EN 10204			Included



### Filling nozzles from the 30-030-xxx series to be applied with gas purging nozzles:

- 0.6 ID up to 3.2 ID filling nozzles can be applied with 30-031-050
- 4.8 ID and 6.0 ID filling nozzles can be used with 30-031-090