

# FP50 SUPPORT

WITH CASTORS



- Designed for support of FP50 - filling and stoppering machine
- Places FP50 in correct working height
- Moveable for saving floor space
- Optional left or right position of control panel
- Option for storage of accessories in lower cabinet

## Watson-Marlow Flexicon

Frejasvej 2-6  
DK-4100 Ringsted  
Denmark  
Tel. +45 57 67 11 55  
Fax. +45 57 67 05 41  
flexicon@flexicon.dk  
www.wmflexicon.dk



# FP50 support TECHNICAL SPECIFICATIONS

The weight and height of FP50 – filling and plugging unit, creates a demand for a lower and stronger support, than offered by an ordinary table.

This support is designed according to the requirements above. Furthermore this support is designed with casters, which allows to move the FP50, when not in use or during cleaning.

The support is made from materials compatible with cleanroom production.

The support offers the choice of whether to place the control panel to the left or the right of the FP50.

If the FP50 is with inlet tray support in straight-line layout, the elevated platform for the control panel can be removed, and the control panel placed below the inlet tray support.

The lower cabinet of the support is designed as a cupboard, for storage of accessories, tools etc.

### Materials and finish:

All stainless steel AISI304

Top plate is polished – other parts are glass blasted

### Weight:

90 kg

