

BioTube Applicator™ Connect with confidence







The BioTube Applicator tube to hose barb assembly machine allows effortless insertion of connectors into a wide variety of tubing, including silicone and thermoplastic elastomers. Union of the tube and connectors is completed without damage to the internal bore of the tube, and importantly, without lubrication such as Isopropyl alcohol (IPA). Providing a repeatable process, this aids cGMP compliance and validation procedures. The presence of IPA in the fluid path is a growing concern in bioprocessing, so avoiding its use simplifies validation.

• **Increase quality and productivity** Automating your tube to hose barb assembly provides repeatability not available with manual processes, thus increasing quality and reducing the risk of junction failures, especially in higher pressure downstream bioprocessing. Automation can also boost productivity and streamline your workflow.

• **Improve safety** Advanced features that can replicate a manual assembly process ensure that even the most complex application can be transferred to the BioTube Applicator. Improving safety in the assembly process, the BioTube Applicator removes risk of operator error and health risks associated with RSI.

• Add flexibility Utilising different combinations of change parts, a wide range of applications can be achieved, ensuring that the needs of your process assembly operation are supported.

The BioTube Applicator utilises change parts to allow different combinations of products to be assembled. These consist of these three components below:

Jaws - patented action ensures the

Tube support - to hold the tubing during the application process

Dio TUBO .





The change parts can be interchanged within 30 seconds and do not require tooling. By utilising change parts we guarantee a repeatable process to ensure every connection is the same as the

Application videos are available which demonstrate the operation of the BioTube Applicator and its use to solve a variety of applications. www.wmfts.com/bta

Cradle - supports the component during the application process



By utilising change parts we guarantee a repeatable process to ensure every connection is the same as the last. This provides a base for validation and aids compliance to cGMP.

Specification







1	Automatic shut down safety cover	
2	Tube support	
3	Jaws	
4	Cradle	
5	Two button activation	
6	Ergonomical, lightweight design including carrying handles	
7	Pressure gauge	
8	Air vent exhaust: fitted with exhaust filter that meets the requirements for use within ISO Class 4 clean room	
9	Air inlet	
10	Power inlet	
11	Pecking dial: can be adjusted to ensure a precise replication of the manual assembly process is achieved	
Model number		BTA125/625
Tubing/component size compatibility		3.2mm - 25.4mm (1/8" - 1")
Main feed, electrical supply		100 - 120V / 220-240V AC 50/60Hz
Control, electrical supply		24V DC
Fuse rating		5 Amps
Compressed air pressure		5- 6 bar
Working temperature		-40C - 70C
Storage temperature		5C - 40C
Weight		18.5kg
Dimensions		220mm(H) x 440mm(W) x 540mm(L)



BIOTECHNOLOGY AND PHARMACEUTICAL SOLUTIONS













Watson-Marlow Fluid Technology Solutions

Watson-Marlow Fluid Technology Solutions supports its customers locally through an extensive global network of direct sales operations and distributors

wmfts.com/global





Disclaimer: The information contained in this document is believed to be correct but Watson-Marlow Limited accepts no liability for any errors it contains and reserves the right to alter specifications without notice. It is the users responsibility to ensure product suitability for use within their application. Watson-Marlow, BioClamp and puresu are registered trademarks of Watson-Marlow Limited. Bio Y, BioBarb, FlatBioEndCap, BioChadcap, BioValve and BioTube applicator are Trademarks of Watson-Marlow Limited. Tri-Clamp is a registered trademark of Ladish Company.

A Spirax-Sarco Engineering plc company