

Accusil Elements

Accusil tubing

Features and benefits

- Excellent dispensing accuracy with $\pm 0.5\%$
- 100 % traceability with laser-etched part-number, lot number and use-by date
- Eight different sizes for filling between 0.1 ml to more than 250 ml

Technical specifications

	Accusil Elements
Sterilisation methods	Autoclave, EtO, Gamma, X-ray
Autoclave stability	121 °C 30 minutes
Gamma stability	Max 50 kGy
Bore size	1.6 - 8 mm
Bore size	0.0625 - 0.315 in
Wall thickness	2.4 mm
Wall thickness	0.09375 in
Colour	Colourless
Transparency	Translucent
Spallation	Low
Operating temperature range	-20 °C to 80 °C
Operating temperature range	-4 °F to 176 °F
Hardness, shore A (5sec)	58 to 64
Specific gravity	1.15 g/ml
Compression set	$\leq 30\%$
Shelf life	5 years

Materials of construction

	Accusil Elements
Tubing	Platinum-cured silicone
Connector material	PVDF

Product Codes

Bore		Wall		Double segment kynar Y-element
mm	inch	mm	inch	
1.6	1/16	2.4	3/32	84-510-016
3.2	1/8	2.4	3/32	84-510-032
4.8	3/16	2.4	3/32	84-510-048
6.4	1/4	2.4	3/32	84-510-064
8.0	5/16	2.4	3/32	84-510-080
9.6	3/8	2.4	3/32	84-510-096

Disclaimer: The information contained in this document is believed to be correct but Watson-Marlow Limited accepts no liability for any errors it contains and reserves the right to alter specifications without notice. It is the user's responsibility to ensure product suitability for use within their application. Watson-Marlow and WMArchitect are registered trademarks of Watson-Marlow Limited. BioClamp, BioBarb, FlatBioEndCap, BioEndCap, BioValve and BioTube applicator are Trademarks of BioPure Technology Limited. Tri-Clamp is a registered trademark of Alfa Laval Corporate AB. Q-Clamp is used under licence.



wmfts.com/global
16 January 2025