Certa Plus 400

Certa Plus process pump



Features and benefits

- Gentle product handling preserves high quality of product
- Impressive suction capability
- Almost zero pulsation
- Outstanding energy efficient pumping principle
- Bi-directional running
- Only one shaft and one seal system
- Different port orientations available, including self-draining position
- Traceability back to base materials












Certa Plus 400 performance

Performance curves

These performance curves show the theoretical flow without consideration of slip




Technical specifications

	Certa Plus 400									
Max. flow rate	41760 L/h									
Max. flow rate	183.2 USGPM									
Volume per revolution	1.16 L									
Volume per revolution	0.305 USG									
Max. particle size	38 mm									
Max. particle size	1.5 in									
Suction pressure	85% of a vacuum									
Max. operating speed	600 rpm									
Max. operating pressure	6 bar									
Max. operating pressure	87 psi									
Operating temperature range	-20 to 100 °C									
Operating temperature range	-4 to 212 °F									
Max. SIP temperature	121 °C									
Max. SIP temperature	250 °F									
Shaft diameter	50 mm									
Shaft diameter	1.97 in									
Shaft height	200 mm									
Shaft height	7.87 in									
Standards	3-A, ATEX, EHEDG Type El Class I, EHEDG Type EL Class I Aseptic									
Certification	EC 1935/2004, FDA, USP ClassVI									
Port orientation	Various options to meet your requirements									
Port type	Flange, Milk fitting, RJT, SMS, TC									

Materials of construction

	Certa Plus 400							
Housing	Stainless steel 1.4435, Stainless steel 316L							
Bearing frame	Stainless steel 1.4301, Stainless steel 304							
Rotor	Stainless steel 1.4435, Stainless steel 316L							
Gate	PEEK, PTFE							
Seals/gaskets	EPDM, FKM							
Shaft seal	PTFE, SSIC							

Certa Plus 400 dimensions



	Nozzles						Foot									Length				Height				Coupling						
Model	1odel N1		N2		N3		F1		F2		F3		F4		F5		L1		L2		H1		H2		C1		C2		С3	
	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.	mm	in.
C400	323.5	12.74	169	6.65	31	1.22	266.5	10.49	254	10.00	Ø21	Ø0.83	301	11.85	169	6.65	513.5	20.22	77	3.03	200	7.87	292.5	11.52	Ø50	Ø1.97	14	0.55	53.5	2.11

Disclaimer: All flow rates shown were obtained pumping water at 20 °C (68 °F) with zero suction and delivery heads. Disclaimer: The information contained in this document is believed to be correct but Watson-Marlow Limited accepts no liability for any errors it contains and reserves the right to alter specifications without notice. Watson-Marlow, MasoSine, Qdos, ReNu, LoadSure, Bioprene, Pumpsil and Marprene are trademarks of Watson-Marlow Limited. STA-PURE PFL® and STA-PURE PCS® are registered trademarks of W.L Gore & Associates Inc. Please state the product code when ordering pumps and tubing.

wmfts.com/global



Monday, May 15, 2023