

# 400/N4 pumphead

400/M and N series pumpheads

## Features and benefits

- 4 channel compact pump design
- Works with 1.0 mm wall thickness tubing and barbed fittings
- Complete pump with motor, or as pumphead only
- Pumphead available in PVDF for greater chemical compatibility




## 400/N4 pumphead performance

Flow rates (ml/min)		
Tube bore (mm)	ml/rev	Max continuous flow (ml/min)
0.5	0.013	7
1.0	0.05	25
2.0	0.18	90
3.0	0.33	150

400/N

## Pressure / suction charts

Santoprene tubing @100rpm



## Technical specifications

	400/N4 pumphead
Number of channels	4
Max. flow rate continuous	150 ml/min
Operating temperature	-10 °C to 45 °C (41 °F to 113 °F)
Pump type	Panel mount
Compatible tubing bore size	0.5, 1, 2, 3 mm
Compatible tubing wall thickness	1 mm

