GORE® STA-PURE Series PCS Elements



GORE® STA-PURE Series PCS

Features and benefits

- Pressure rated up to 7 bar (100 psi)
- Longest available tube life
- STA-PURE PCS's bore is three times smoother than thermoplastic elastomer tubing, with no grooves or cracking developing as flexing continues
- Spallation virtually eliminated, leaving high-value duty fluids uncontaminated
- Available in 20 different bore sizes, three wall thicknesses and three element lengths
- Operating temperature -20 °C to 80 °C
- STA-PURE PCS is a trademark of WL Gore & Associates inc.



Technical specifications

	GORE® STA-PURE Series PCS Elements	
Bore size	1.6 - 25.4 mm	
Bore size	0.0625 - 1 in	
Wall thickness	2.4 - 4.8 mm	
Wall thickness	0.0938 - 0.1875 in	
Colour	Off-white	
Transparency	Opaque	
Spallation	Very low	
Certification	USP <87>, USP <88>, USP Class VI	
Operating temperature range	-20 °C to 80 °C	
Operating temperature range	-4 °F to 176 °F	
Sterilisation methods	Autoclave, CIP, SIP	

Materials of construction

	GORE® STA-PURE Series PCS Elements
Tubing	ePTFE and platinum-cured silicone composite
Connector material	PVDF

Product codes

STA-PURE PCS Double Segment Kynar Y-Element ordering codes						
Bore		Wall Thickn	ess	Product codes		
mm	inch	mm	inch			
1.6	0.0625	2.4	0.0938	961.E016.K24		
3.2	0.1250	2.4	0.0938	961.E032.K24		
4.8	0.1875	2.4	0.0938	961.E048.K24		
6.4	0.2500	2.4	0.0938	961.E064.K24		
8.0	0.3125	2.4	0.0938	961.E080.K24		
8.0	0.3125	4.0	0.1562	961.E080.K40		
12.0	0.4688	2.4	0.0938	961.E120.K40		
16.0	0.6250	2.4	0.0938	961.E160.K40		

STA-PURE PCS Loadsure Elements ordering codes				
Product Name	Product codes			
STA-PURE PCS 520 LoadSure 3.2mm Bore Tri-Clamp Element	961.0032.PFT			
STA-PURE PCS 520 LoadSure 6.4mm Bore Tri-Clamp Element	961.0064.PFT			
STA-PURE PCS 520 LoadSure 9.6mm Bore Tri-Clamp Element	961.0096.PFT			
STA-PURE PCS 620 LoadSure 12.0mm Bore Tri-Clamp Element	961.0120.PFT			
STA-PURE PCS 620 LoadSure 17.0mm Bore Tri-Clamp Element	961.0170.PFT			
STA-PURE PCS 720 LoadSure 12.7mm Bore Tri-Clamp Element	961.0127.PFT			
STA-PURE PCS 720 LoadSure 15.9mm Bore Tri-Clamp Element	961.0159.PFT			
STA-PURE PCS 720 LoadSure 19.0mm Bore Tri-Clamp Element	961.0190.PFT			
STA-PURE PCS 720 LoadSure 25.4mm Bore Tri-Clamp Element	961.0254.PFT -			

Disclaimer: Watson-Marlow and LoadSure are registered trademarks of Watson-Marlow Limited. STA-PURE PCS is trademark of WL Gore & Associates Inc. Please state the product code when ordering pumps and tubing.

wmfts.com/global



22 November 2023