


Fluid Transfer Tubing and hose sets

Single-use fluid transfer assemblies

Features and benefits

- Complete product traceability
- Large component library to choose from
- Support of internal validation protocols and validation guides available upon request
- Standard configurations available for reduced turn-around times
- Standard assemblies sold non-irradiated, irradiated and sterile assemblies on request
- Flexible integration with an open architecture approach
- Use in combination with Watson-Marlow peristaltic pumps for an end-to-end fluid management solution
- Supply security with Watson-Marlow's global manufacturing and supply network




Technical specifications

	Fluid Transfer Tubing and hose sets
Certification and compliance	ADCF, ISO10993, USP <788>, USP <85>, USP Class VI
Manufacturing	ISO 14644-1 Class 7, ISO 9001
Sterility	Non-irradiated (irradiated or sterile on request)

Materials of construction

	Fluid Transfer Tubing and hose sets
Tubing	SEBS, Silicone, Thermoplastic elastomer
Connector material	Polycarbonate (PC), Polypropylene, Polysulfone (UDEL 1700), PVDF

Order Information



Part Number	Part Name (Description)	Pumps	Pumphead
33-0913-211112	Pumpsil Kit - 3.2x1.6mm - .3m long w/TC-NON	500 Series	520R
33-0913-211113	Pumpsil Kit - 3.2x1.6mm - .3m long w/MPC-NON	500 Series	520R
33-0913-211114	Pumpsil Kit - 3.2x1.6mm - .3m long w/Aseptiquik-NON	500 Series	520R
33-0913-211122	Pumpsil Kit - 6.4x1.6mm - .3m long w/TC-NON	Transfer use only	
33-0913-211123	Pumpsil Kit - 6.4x1.6mm - .3m long w/MPC-NON	Transfer use only	
33-0913-211124	Pumpsil Kit - 6.4x1.6mm - .3m long w/Aseptiquik-NON	Transfer use only	
33-0913-211212	Pumpsil Kit - 3.2x1.6mm - 1m long w/TC-NON	500 Series	520R
33-0913-211213	Pumpsil Kit - 3.2x1.6mm - 1m long w/MPC-NON	500 Series	520R
33-0913-211214	Pumpsil Kit - 3.2x1.6mm - 1m long w/Aseptiquik-NON	500 Series	520R
33-0913-211222	Pumpsil Kit - 6.4x1.6mm - 1m long w/TC-NON	500 Series	520R
33-0913-211223	Pumpsil Kit - 6.4x1.6mm - 1m long w/MPC-NON	500 Series	520R
33-0913-211224	Pumpsil Kit - 6.4x1.6mm - 1m long w/Aseptiquik-NON	500 Series	520R

33-0913-211232	Pumpsil Kit - 9.6x3.2mm - 1m long w/TC-NON	600 Series	620R
33-0913-211233	Pumpsil Kit - 9.6x3.2mm - 1m long w/MPC-NON	600 Series	620R
33-0913-211234	Pumpsil Kit - 9.6x3.2mm - 1m long w/Aseptiquik-NON	600 Series	620R
33-0913-211242	Pumpsil Kit - 12.7x3.2mm - 1m long w/TC-NON	600 Series	620R
33-0913-211243	Pumpsil Kit - 12.7x3.2mm - 1m long w/MPC-NON	600 Series	620R
33-0913-211244	Pumpsil Kit - 12.7x3.2mm - 1m long w/Aseptiquik-NON	600 Series	620R
33-0913-211252	Pumpsil Kit - 19x4.8mm - 1m long w/TC-NON	700 Series	720R
33-0913-211253	Pumpsil Kit - 19x4.8mm - 1m long w/MPC-NON	700 Series	720R
33-0913-211254	Pumpsil Kit - 19x4.8mm - 1m long w/Aseptiquik-NON	700 Series	720R
33-0913-211312	Pumpsil Kit - 3.2x1.6mm - 2.5m long w/TC-NON	500 Series	520R
33-0913-211313	Pumpsil Kit - 3.2x1.6mm - 2.5m long w/MPC-NON	500 Series	520R
33-0913-211314	Pumpsil Kit - 3.2x1.6mm - 2.5m long w/Aseptiquik-NON	500 Series	520R
33-0913-211322	Pumpsil Kit - 6.4x1.6mm - 2.5m long w/TC-NON	500 Series	520R
33-0913-211323	Pumpsil Kit - 6.4x1.6mm - 2.5m long w/MPC-NON	500 Series	520R
33-0913-211324	Pumpsil Kit - 6.4x1.6mm - 2.5m long w/Aseptiquik-NON	500 Series	520R
33-0913-211332	Pumpsil Kit - 9.6x3.2mm - 2.5m long w/TC-NON	600 Series	620R
33-0913-211333	Pumpsil Kit - 9.6x3.2mm - 2.5m long w/MPC-NON	600 Series	620R
33-0913-211334	Pumpsil Kit - 9.6x3.2mm - 2.5m long w/Aseptiquik-NON	600 Series	620R
33-0913-211342	Pumpsil Kit - 12.7x3.2mm - 2.5m long w/TC-NON	600 Series	620R
33-0913-211343	Pumpsil Kit - 12.7x3.2mm - 2.5m long w/MPC-NON	600 Series	620R
33-0913-211344	Pumpsil Kit - 12.7x3.2mm - 2.5m long w/Aseptiquik-NON	600 Series	620R
33-0913-211352	Pumpsil Kit - 19x4.8mm - 2.5m long w/TC-NON	700 Series	720R
33-0913-211353	Pumpsil Kit - 19x4.8mm - 2.5m long w/MPC-NON	700 Series	720R
33-0913-211354	Pumpsil Kit - 19x4.8mm - 2.5m long w/Aseptiquik-NON	700 Series	720R
33-0913-211432	Pumpsil Kit - 9.6x3.2mm - 5m long w/TC-NON	600 Series	620R
33-0913-211433	Pumpsil Kit - 9.6x3.2mm - 5m long w/MPC-NON	600 Series	620R
33-0913-211434	Pumpsil Kit - 9.6x3.2mm - 5m long w/Aseptiquik-NON	600 Series	620R
33-0913-211442	Pumpsil Kit - 12.7x3.2mm - 5m long w/TC-NON	600 Series	620R
33-0913-211443	Pumpsil Kit - 12.7x3.2mm - 5m long w/MPC-NON	600 Series	620R
33-0913-211444	Pumpsil Kit - 12.7x3.2mm - 5m long w/Aseptiquik-NON	600 Series	620R
33-0913-211452	Pumpsil Kit - 19x4.8mm - 5m long w/TC-NON	700 Series	720R
33-0913-211453	Pumpsil Kit - 19x4.8mm - 5m long w/MPC-NON	700 Series	720R
33-0913-211454	Pumpsil Kit - 19x4.8mm - 5m long w/Aseptiquik-NON	700 Series	720R
33-0933-211113	Bioprene Kit - 3.2x1.6mm - .3m long w/MPC-NON	Transfer use only	
33-0933-211114	Bioprene Kit - 3.2x1.6mm - .3m long w/Aseptiquik-NON	Transfer use only	
33-0933-211122	Bioprene Kit - 6.4x1.6mm - .3m long w/TC-NON	Transfer use only	
33-0933-211123	Bioprene Kit - 6.4x1.6mm - .3m long w/MPC-NON	Transfer use only	
33-0933-211124	Bioprene Kit - 6.4x1.6mm - .3m long w/Aseptiquik-NON	Transfer use only	
33-0933-211213	Bioprene Kit - 3.2x1.6mm - 1m long w/MPC-NON	500 Series	520R
33-0933-211214	Bioprene Kit - 3.2x1.6mm - 1m long w/Aseptiquik-NON	500 Series	520R
33-0933-211222	Bioprene Kit - 6.4x1.6mm - 1m long w/TC-NON	500 Series	520R
33-0933-211223	Bioprene Kit - 6.4x1.6mm - 1m long w/MPC-NON	500 Series	520R
33-0933-211224	Bioprene Kit - 6.4x1.6mm - 1m long w/Aseptiquik-NON	500 Series	520R
33-0933-211232	Bioprene Kit - 9.6x3.2mm - 1m long w/TC-NON	600 Series	620R
33-0933-211233	Bioprene Kit - 9.6x3.2mm - 1m long w/MPC-NON	600 Series	620R
33-0933-211234	Bioprene Kit - 9.6x3.2mm - 1m long w/Aseptiquik-NON	600 Series	620R
33-0933-211242	Bioprene Kit - 12.7x3.2mm - 1m long w/TC-NON	600 Series	620R
33-0933-211243	Bioprene Kit - 12.7x3.2mm - 1m long w/MPC-NON	600 Series	620R
33-0933-211244	Bioprene Kit - 12.7x3.2mm - 1m long w/Aseptiquik-NON	600 Series	620R
33-0933-211252	Bioprene Kit - 19x4.8mm - 1m long w/TC-NON	700 Series	720R
33-0933-211253	Bioprene Kit - 19x4.8mm - 1m long w/MPC-NON	700 Series	720R
33-0933-211254	Bioprene Kit - 19x4.8mm - 1m long w/Aseptiquik-NON	700 Series	720R
33-0933-211313	Bioprene Kit - 3.2x1.6mm - 2.5m long w/MPC-NON	500 Series	520R
33-0933-211314	Bioprene Kit - 3.2x1.6mm - 2.5m long w/Aseptiquik-NON	500 Series	520R
33-0933-211322	Bioprene Kit - 6.4x1.6mm - 2.5m long w/TC-NON	500 Series	520R
33-0933-211323	Bioprene Kit - 6.4x1.6mm - 2.5m long w/MPC-NON	500 Series	520R
33-0933-211324	Bioprene Kit - 6.4x1.6mm - 2.5m long w/Aseptiquik-NON	500 Series	520R
33-0933-211332	Bioprene Kit - 9.6x3.2mm - 2.5m long w/TC-NON	600 Series	620R
33-0933-211333	Bioprene Kit - 9.6x3.2mm - 2.5m long w/MPC-NON	600 Series	620R
33-0933-211334	Bioprene Kit - 9.6x3.2mm - 2.5m long w/Aseptiquik-NON	600 Series	620R
33-0933-211342	Bioprene Kit - 12.7x3.2mm - 2.5m long w/TC-NON	600 Series	620R
33-0933-211343	Bioprene Kit - 12.7x3.2mm - 2.5m long w/MPC-NON	600 Series	620R
33-0933-211344	Bioprene Kit - 12.7x3.2mm - 2.5m long w/Aseptiquik-NON	600 Series	620R
33-0933-211352	Bioprene Kit - 19x4.8mm - 2.5m long w/TC-NON	700 Series	720R
33-0933-211353	Bioprene Kit - 19x4.8mm - 2.5m long w/MPC-NON	700 Series	720R
33-0933-211354	Bioprene Kit - 19x4.8mm - 2.5m long w/Aseptiquik-NON	700 Series	720R
33-0933-211432	Bioprene Kit - 9.6x3.2mm - 5m long w/TC-NON	600 Series	620R
33-0933-211433	Bioprene Kit - 9.6x3.2mm - 5m long w/MPC-NON	600 Series	620R
33-0933-211434	Bioprene Kit - 9.6x3.2mm - 5m long w/Aseptiquik-NON	600 Series	620R
33-0933-211442	Bioprene Kit - 12.7x3.2mm - 5m long w/TC-NON	600 Series	620R
33-0933-211443	Bioprene Kit - 12.7x3.2mm - 5m long w/MPC-NON	600 Series	620R
33-0933-211444	Bioprene Kit - 12.7x3.2mm - 5m long w/Aseptiquik-NON	600 Series	620R
33-0933-211452	Bioprene Kit - 19x4.8mm - 5m long w/TC-NON	700 Series	720R
33-0933-211453	Bioprene Kit - 19x4.8mm - 5m long w/MPC-NON	700 Series	720R
33-0933-211454	Bioprene Kit - 19x4.8mm - 5m long w/Aseptiquik-NON	700 Series	720R
33-0942-211112	Pureweld Kit - 3.2x1.6mm - .3m long w/TC-NON	Transfer use only	
33-0942-211113	Pureweld Kit - 3.2x1.6mm - .3m long w/MPC-NON	Transfer use only	
33-0942-211114	Pureweld Kit - 3.2x1.6mm - .3m long w/Aseptiquik-NON	Transfer use only	
33-0942-211122	Pureweld Kit - 6.4x1.6mm - .3m long w/TC-NON	Transfer use only	
33-0942-211123	Pureweld Kit - 6.4x1.6mm - .3m long w/MPC-NON	Transfer use only	
33-0942-211124	Pureweld Kit - 6.4x1.6mm - .3m long w/Aseptiquik-NON	Transfer use only	
33-0942-211212	Pureweld Kit - 3.2x1.6mm - 1m long w/TC-NON	500 Series	520R
33-0942-211213	Pureweld Kit - 3.2x1.6mm - 1m long w/MPC-NON	500 Series	520R
33-0942-211214	Pureweld Kit - 3.2x1.6mm - 1m long w/Aseptiquik-NON	500 Series	520R
33-0942-211222	Pureweld Kit - 6.4x1.6mm - 1m long w/TC-NON	500 Series	520R
33-0942-211223	Pureweld Kit - 6.4x1.6mm - 1m long w/MPC-NON	500 Series	520R
33-0942-211224	Pureweld Kit - 6.4x1.6mm - 1m long w/Aseptiquik-NON	500 Series	520R
33-0942-211232	Pureweld Kit - 9.6x3.2mm - 1m long w/TC-NON	600 Series	620R

33-0942-211233	Pureweld Kit - 9.6x3.2mm - 1m long w/MPC-NON	600 Series	620R
33-0942-211234	Pureweld Kit - 9.6x3.2mm - 1m long w/Aseptiquik-NON	600 Series	620R
33-0942-211242	Pureweld Kit - 12.7x3.2mm - 1m long w/TC-NON	600 Series	620R
33-0942-211243	Pureweld Kit - 12.7x3.2mm - 1m long w/MPC-NON	600 Series	620R
33-0942-211244	Pureweld Kit - 12.7x3.2mm - 1m long w/Aseptiquik-NON	600 Series	620R
33-0942-211252	Pureweld Kit - 19x4.8mm - 1m long w/TC-NON	700 Series	720R
33-0942-211253	Pureweld Kit - 19x4.8mm - 1m long w/MPC-NON	700 Series	720R
33-0942-211254	Pureweld Kit - 19x4.8mm - 1m long w/Aseptiquik-NON	700 Series	720R
33-0942-211312	Pureweld Kit - 3.2x1.6mm - 2.5m long w/TC-NON	500 Series	520R
33-0942-211313	Pureweld Kit - 3.2x1.6mm - 2.5m long w/MPC-NON	500 Series	520R
33-0942-211314	Pureweld Kit - 3.2x1.6mm - 2.5m long w/Aseptiquik-NON	500 Series	520R
33-0942-211322	Pureweld Kit - 6.4x1.6mm - 2.5m long w/TC-NON	500 Series	520R
33-0942-211323	Pureweld Kit - 6.4x1.6mm - 2.5m long w/MPC-NON	500 Series	520R
33-0942-211324	Pureweld Kit - 6.4x1.6mm - 2.5m long w/Aseptiquik-NON	500 Series	520R
33-0942-211332	Pureweld Kit - 9.6x3.2mm - 2.5m long w/TC-NON	600 Series	620R
33-0942-211333	Pureweld Kit - 9.6x3.2mm - 2.5m long w/MPC-NON	600 Series	620R
33-0942-211334	Pureweld Kit - 9.6x3.2mm - 2.5m long w/Aseptiquik-NON	600 Series	620R
33-0942-211342	Pureweld Kit - 12.7x3.2mm - 2.5m long w/TC-NON	600 Series	620R
33-0942-211343	Pureweld Kit - 12.7x3.2mm - 2.5m long w/MPC-NON	600 Series	620R
33-0942-211344	Pureweld Kit - 12.7x3.2mm - 2.5m long w/Aseptiquik-NON	600 Series	620R
33-0942-211352	Pureweld Kit - 19x4.8mm - 2.5m long w/TC-NON	700 Series	720R
33-0942-211353	Pureweld Kit - 19x4.8mm - 2.5m long w/MPC-NON	700 Series	720R
33-0942-211354	Pureweld Kit - 19x4.8mm - 2.5m long w/Aseptiquik-NON	700 Series	720R
33-0942-211432	Pureweld Kit - 9.6x3.2mm - 5m long w/TC-NON	600 Series	620R
33-0942-211433	Pureweld Kit - 9.6x3.2mm - 5m long w/MPC-NON	600 Series	620R
33-0942-211434	Pureweld Kit - 9.6x3.2mm - 5m long w/Aseptiquik-NON	600 Series	620R
33-0942-211442	Pureweld Kit - 12.7x3.2mm - 5m long w/TC-NON	600 Series	620R
33-0942-211443	Pureweld Kit - 12.7x3.2mm - 5m long w/MPC-NON	600 Series	620R
33-0942-211444	Pureweld Kit - 12.7x3.2mm - 5m long w/Aseptiquik-NON	600 Series	620R
33-0942-211452	Pureweld Kit - 19x4.8mm - 5m long w/TC-NON	700 Series	720R
33-0942-211453	Pureweld Kit - 19x4.8mm - 5m long w/MPC-NON	700 Series	720R
33-0942-211454	Pureweld Kit - 19x4.8mm - 5m long w/Aseptiquik-NON	700 Series	720R
33-BSHP-211214	Braided Hose Kit-3/8"ID-1m long-TC-NON	Transfer use only	
33-BSHP-211222	Braided Hose Kit-3/8"ID-1m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211224	Braided Hose Kit-1/2"ID-1m long-TC-NON	Transfer use only	
33-BSHP-211232	Braided Hose Kit-1/2"ID-1m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211234	Braided Hose Kit-3/4"ID-1m long-TC-NON	Transfer use only	
33-BSHP-211242	Braided Hose Kit-3/4"ID-1m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211244	Braided Hose Kit-1"ID-1m long-TC-NON	Transfer use only	
33-BSHP-211312	Braided Hose Kit-1"ID-1m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211314	Braided Hose Kit-3/8"ID-2.5m long-TC-NON	Transfer use only	
33-BSHP-211322	Braided Hose Kit-3/8"ID-2.5m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211324	Braided Hose Kit-1/2"ID-2.5m long-TC-NON	Transfer use only	
33-BSHP-211332	Braided Hose Kit-1/2"ID-2.5m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211334	Braided Hose Kit-3/4"ID-2.5m long-TC-NON	Transfer use only	
33-BSHP-211342	Braided Hose Kit-3/4"ID-2.5m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211344	Braided Hose Kit-1"ID-2.5m long-TC-NON	Transfer use only	
33-BSHP-211412	Braided Hose Kit-1"ID-2.5m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211414	Braided Hose Kit-3/8"ID-5m long-TC-NON	Transfer use only	
33-BSHP-211422	Braided Hose Kit-3/8"ID-5m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211424	Braided Hose Kit-1/2"ID-5m long-TC-NON	Transfer use only	
33-BSHP-211432	Braided Hose Kit-1/2"ID-5m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211434	Braided Hose Kit-3/4"ID-5m long-TC-NON	Transfer use only	
33-BSHP-211442	Braided Hose Kit-3/4"ID-5m long-Aseptiquik-NON	Transfer use only	
33-BSHP-211444	Braided Hose Kit-1"ID-5m long-TC-NON	Transfer use only	

Standard configurations have been designed based on the most used components. Assemblies listed are non-irradiated. Triclamp fittings are uncapped and standard size. Detailed drawings with bill of materials available upon request. For any modifications to length, fitting size/type or irradiation requirement please contact your local Watson-Marlow representative.

Disclaimer: The information contained in this document is believed to be correct but Watson-Marlow Limited accepts no liability for any errors it contains and reserves the right to alter specifications without notice. It is the user's responsibility to ensure product suitability for use within their application. Watson-Marlow and WMArchitect are registered trademarks of Watson-Marlow Limited. BioClamp, BioBarb, FlatBioEndCap, BioEndCap, BioValve, BioTube applicator and Q-Clamp are Trademarks of BioPure Technology Limited. Tri-Clamp is a registered trademark of Alfa Laval Corporate AB.



wmfms.com/global
01 August 2025